SDC-G5-UltraSoft Instructions for use – ENG

Items reference:

- 1106 (G5-UltraSoft 0.09mm, two sides diamond coated)

To be used in combination with the following connection accessories, items reference:

- 1990 (G5-Shank pink)
- FS-150 (G5-Shank Ohead blue)

Compatible with all EVA oscillating system hand-pieces on the market with an instrument insertion diameter $\varnothing=2.30\div 2.35$ mm.

Particularly indicated to be used with the following hand-pieces, items reference:

- 1971 (NSK Ti-Max X55L, optic)
- 1972 (NSK Ti-Max X55, non-optic)
- 3001 (Nose Cone) + 1920 (SDC O-Head)

1. Purpose

High gloss polishing of the enamel surface, composites and/or cements after inter-proximal reduction, restorations and/or prophylaxis treatments.

2. Introduction

For traditional dentists and orthodontics is not usual to polish enamel of the treated area after inter-proximal reduction with a fine grit file of $6\mu\text{m}$. Normally inter-proximal reduction treatments were finished with a $15\mu\text{m}$ grit file.

Maintaining the treated surface as smooth as the natural one is fundamental for the oral cavity health.

3. General information

Oscillating instrument, consisting of a calibrated file, used for a controlled high gloss polishing of enamel surface, composites and/or cements after inter-proximal reduction, restorations and/or prophylaxis treatments.

- Oscillating diamond coated metal strip
- Diamond grit sizes: $6\mu\text{m}$
- Thickness of the strips: 0.09mm
- Height of the strip: 4mm
- Length of the working part: 18mm
- Sterilisable and re-usable

The guidelines implemented and regularly applied for the design and production of this file device, are imposed by the following international standards:

- UNI CEI EN 1041:2009
- UNI EN ISO 1639:2010
- UNI CEI EN ISO 15223-1:2012
- UNI CEI EN ISO 14971:2012
- UNI EN ISO 10993-1:2010
- IEC 62366-1:2015

4. How to Install connection accessory and file

- Insert the shank (connection accessory) from its end into the front part of the hand-piece head taking care to choose the right shank indifferently of the supplied hand-piece:
 - All EVA oscillating system hand-pieces: pink shank (item ref. 1990)
 - SDC O-Head: blue shank (item ref. FS-150)
- Slide the shank until it is fully inserted and verify it is locked into the hand-piece
- At the bottom side of the shank there is a little stopper which prevents the instrument from slipping out: insert the file from the bottom side into the shank and push upwards until it clicks fit

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5. Operating Principles

- Once the tool is inserted into the hand-piece head be sure it is fully and firmly fixed
- Switch on the hand-piece by turning on the micro-motor
- Set and bring the instrument to the recommended speed before applying it on the tooth
- With the file in motion enter the inter-proximal space to be treated with a cervical motion
- Proceed with the treatment paying attention to the speed and the pressure applied
- Once the treatment is done remove the instrument still in movement from interdental space
- Switch off the oscillating hand-piece by turning off the micro-motor

6. Principles of safety

- Provide sufficient irrigation: at least 50 mL / min
- Do not exceed the recommended speed
 - Hand-piece NSK Ti-Max X55 / X55L: 40'000 rpm = 10'000 stk. / min,
 - Red Nose Cone + SDC O-Head: 20'000 rpm = 5'000 stk. / min
 - All other EVA oscillating system hand-pieces: refer to manufacturer instructions
- Do not exceed the recommended pressure: 1N ÷ 3N (max. 3N)
- Do not stop the motion of the file while the same is in the interdental space
- In case of overcrowding open the interdental space using a wedge in the cervical area
- Do not bend, warp or leverage the instrument
- Not observing these precautions leads to a high risk in terms of safety and can lead to the possible breakage or damage of the instrument
- SDC accepts no liability for improper use of its products

7. Warnings and Precautions

- The instrument is for exclusive professional use, mainly for dentists specialized in Orthodontics 
- Do not bend, warp or leverage the instrument over the limit: the use of excessive force (leverage with the hand-piece) can lead to premature breakage
- Excessive lateral swerve or obstruction during the oscillation movement increases the risk of instrument damage
- Do not exceed the 20'000 oscillations per minute to avoid breakage of the instrument
- Temperatures over 180°C are to be avoided in order to not reduce the instrument lifetime
- Avoid forces exceeding 3N because in this case you risk damaging the instrument
- The consumption of the instrument lifetime until its worn out, replacement and disposal is an exclusive decision of the physician according to his operative needs
- Verify instrument integrity before reutilization: absence of small cracks, breakages and/or deformations on plastic part, diamond presence on metal band
- Be careful not to injure the soft tissues, especially when you cross the equator of the tooth
- To ensure traceability of the instrument during the entire instrument lifetime we recommend to preserve the packaging
- On the product packages SDC lists the references that allow you to trace the entire production cycle of each product. In case of detected problems of a constructive nature of SDC product, we ask you to return in its full packaging, in order to allow us to be able to trace the cause of the problem that we will report

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8. Maintenance

- The instrument is supplied non-sterile: sterilize before the first use
- In order to preserve the abrasive power of the instrument, after each use the instrument must be cleaned of any residues, disinfected and sterilized

9. Cleaning

- Disinfection can take place in an ultrasonic cleaning device containing a liquid degreaser and / or a chemically inhibited disinfectant
- The cleaning of the instrument after each use can be done by brushing the diamond part provided that the brush bristles are not too rigid to avoid damaging metal strip, diamond and polymer

10. Sterilization

- Use single or fractionated vacuum autoclave and vacuum drying
- The hot air sterilizers are not suitable for diamond instruments
- Sterilization must follow validated procedures in compliance with European standard EN 13060:
 - B Standard 134 (T=135.5°C, P=2.16bar, t_{Plateau}=4', t_{Drying}=15', t_{Total}=30'÷40')
 - B Prion 134 (T=135.5°C, P=2.16bar, t_{Plateau}=18', t_{Drying}=15', t_{Total}=44'÷54')
 - B Standard 121 (T=122.5°C, P=1.14bar, t_{Plateau}=15', t_{Drying}=20', t_{Total}=50'÷60')

11. Contraindications

- Use only as described above to avoid incurring drawbacks due to the carelessness and distraction
- Do not use on patients suffering from indisposition and / or diseases of any nature and / or type
- Do not disperse in the environment
- The instrument gets in contact with biological tissues during its usage: dispose of the instrument in accordance with the destination country legislation
- Keep out of the sight and reach of children
- Do not swallow

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