SDC-Measuring Gauge Instructions for use – ENG

Items references:

- 1991 (Measuring Gauge with 5 thicknesses)
- 1992 (Measuring Gauge with 9 thicknesses)

1. Purpose

The Swissdentacare Measuring Gauge is an interdental space control device during IPR treatments (Inter Proximal Reduction). It is a manual tool consisting of different thicknesses of stainless steel with ergonomic grip and stainless-steel assembly/disassembly screw to facilitate sterilization operations.

2. Introduction

In orthodontics, when using non-calibrated files for the reduction of interproximal enamel, it is necessary to check how much space has been created between a tooth and its adjacent one. The clinician can determine with the gauge whether the existing or created space is insufficient, adequate, or excessive compared to the guidelines established by existing IPR (Inter Proximal Reduction) technique protocols.

3. General information

The Swissdentacare Measuring Gauges have the following dimensional characteristics:

- 1991 (Measuring Gauge with 5 thicknesses)
 - Thicknesses: 0.10mm, 0.20mm, 0.30mm, 0.40mm, 0.50mm
 - Height of measuring section: 5.50mm
 - Grip height: 17mm
 - Length of measuring section: 23mm
 - Total length of instrument: 48.5mm
 - Sterilisable and re-usable
- 1992 (Measuring Gauge with 9 thicknesses)
 - Thicknesses: 0.10mm, 0.12mm, 0.15mm, 0.20mm, 0.25mm, 0.30mm, 0.35mm, 0.40mm, 0.50mm
 - Height of measuring section: 5.50mm
 - Grip height: 17mm
 - Length of measuring section: 23mm
 - Total length of instrument: 48.5mm
 - Sterilisable and re-usable

The guidelines implemented and regularly applied for the design and production of this instrument, are imposed by the harmonized international standards.

4. How to Install the instrument

Unscrew the fixing screw, clean, disinfect and sterilize all parts of the instrument before first use and after each patient. Reassemble the gauge for next use.

5. Operating Principles

Enter the gauge thickness of the size corresponding to the space created and check if it fits into the existing gap:

- If the thickness does not enter into the existing gap, the space created is insufficient.
- If the shim enters the slot with slight interference, the space created is the correct size.
- Excessive space has been created if the shim enters easily

6. Principles of safety

- Avoid bending, warping, and/or levering with the instrument
- Check the condition of each component of the gauge before assembly and use for damages of the gauge smooth surface, traces, residues, rust, or folds from affecting its safe use.
- Failing to observe these precautions leads to a high risk in terms of safety and can lead to the eventual instrument breakage.

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- Swissdentacare SA accepts no liability for improper use of its products

7. Warnings and Precautions

- The instrument is for the exclusive use of dental professionals, mainly dentists specializing in Orthodontics 
- Do not bend, warp, or leverage the instrument beyond the limit: the use of excessive force (leverage effect) can lead to premature breakage
- Temperatures over 180°C are to be avoided in order to not reduce the instrument lifetime
- Verify instrument integrity before reutilization: absence of small cracks, breakages and/or deformations on metal part, rust.
- Be careful not to injure the soft tissues, especially when you cross the equator of the tooth

8. Maintenance

- The instrument is supplied non-sterile: sterilize before the first use
- To preserve the smooth surface of the gauges, after each use the instrument must be cleaned of any residues, disinfected, and sterilized

9. Cleaning

- Disinfection can take place in an ultrasonic cleaning device containing a liquid degreaser and / or a chemically inhibited disinfectant
- The cleaning of the instrument after each use can be done by using a damp cloth.

10. Sterilization

- Use single or fractionated vacuum autoclave and vacuum drying
- The hot air sterilizers are not suitable
- Sterilization must follow validated procedures in compliance with European standard EN 13060:
 - B Standard 134 (T=135.5°C, P=2.16bar, t_{Plateau}=4', t_{Drying}=15', t_{Total}=30'÷40')
 - B Prion 134 (T=135.5°C, P=2.16bar, t_{Plateau}=18', t_{Drying}=15', t_{Total}=44'÷54')
 - B Standard 121 (T=122.5°C, P=1.14bar, t_{Plateau}=15', t_{Drying}=20', t_{Total}=50'÷60')

11. Contraindications

- Use only as described above to avoid incurring drawbacks due to the carelessness and distraction
- Do not use on patients suffering from indisposition and / or diseases of any nature and / or type
- Do not disperse in the environment
- Keep out of the sight and reach of children
- Do not swallow

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