SDC-Proxocare Instructions for use – ENG

Items reference:

- 1706 (SDC-Proxocare 6µm, one side diamond coated)
- 1715 (SDC-Proxocare 15µm, one side diamond coated)
- 1725 (SDC-Proxocare 25µm, one side diamond coated)
- 1740 (SDC-Proxocare 40µm, one side diamond coated)
- 1760 (SDC-Proxocare 60µm, one side diamond coated)
- 1790 (SDC-Proxocare 90µm, one side diamond coated)

Compatible with all EVA oscillating system hand-pieces on the market with an instrument insertion diameter $\varnothing=2.30\div 2.35$ mm.

Particularly indicated to be used with the following hand-pieces, items reference:

- 1971 (NSK Ti-Max X55L, optic)
- 1972 (NSK Ti-Max X55, non-optic)

1. Purpose

Oscillating diamond coated files particularly suitable for fine preparations and polishing of proximal and concave restorations. Removal of overhanging sections of composites and cements.

2. Introduction

In traditional dentistry it is common to use rotary instruments which are unsafe, inaccurate and very invasive in order to prepare and polish restorations.

The benefits of using this product family are:

- Easier access to the inter-proximal space as the file is thinner
- No accidental damage to adjacent teeth (one side diamond coated file)
- Avoiding generating waves on the threatened site (e.g. caused by rotating instruments)

3. General information

The oscillating files instrument family are used for fine preparations and polishing of restoration treatments.

- Oscillating diamond coated metal file
- Diamond grit sizes: 6µm, 15µm, 25µm, 40µm, 60µm
- Length of the file: 8mm
- Sterilisable and re-usable

The guidelines implemented and regularly applied for the design and production of these file devices, are imposed by the following international standards:

- UNI CEI EN 1041:2009
- UNI EN ISO 1639:2010
- UNI CEI EN ISO 15223-1:2012
- UNI CEI EN ISO 14971:2012
- UNI EN ISO 10993-1:2010
- IEC 62366-1:2015

4. How to Install instrument

- Insert the instrument by the end of its shank into the front part of the hand-piece head
- Push all the way the instrument until it is fully inserted and verify it is locked into the hand-piece head

5. Operating Principles

- Once the tool is inserted into the hand-piece head be sure it is fully and firmly fixed
- Switch on the hand-piece by turning on the micro-motor
- Set and bring the instrument to the recommended speed before applying it on the tooth
- With the instrument in motion enter the inter-proximal space to be treated with a cervical motion
- Proceed with the treatment paying attention to the speed and the pressure applied
- Once the treatment is done remove the instrument still in movement from interdental space
- Switch off the oscillating hand-piece by turning off the micro-motor

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6. Principles of safety

- Provide sufficient irrigation: at least 50 mL / min
- Do not exceed the recommended speed
 - Hand-piece NSK Ti-Max X55 / X55L: 40'000 rpm = 10'000 stk. / min,
 - All other EVA oscillating system hand-pieces: refer to manufacturer instructions
- Do not exceed the recommended pressure: 1N ÷ 3N (max. 3N)
- Do not stop the motion of the instrument while the same is in the interdental space
- In case of overcrowding open the interdental space using a wedge in the cervical area
- Do not bend, warp or leverage the instrument
- Not observing these precautions leads to a high risk in terms of safety and can lead to the possible breakage or damage of the instrument
- SDC accepts no liability for improper use of its products

7. Warnings and Precautions

- The instrument is for exclusive professional use, mainly for dentists specialized in Restorative Dentistry 
- Do not bend, warp or leverage the instrument over the limit: the use of excessive force (leverage with the hand-piece) can lead to premature breakage
- Excessive lateral swerve or obstruction during the oscillation movement increases the risk of instrument damage
- Do not exceed the 20'000 oscillations per minute to avoid breakage of the instrument
- Temperatures over 180°C are to be avoided in order to not reduce the instrument lifetime
- Avoid forces exceeding 3N because in this case you risk damaging the instrument
- The consumption of the instrument lifetime until its worn out, replacement and disposal is an exclusive decision of the physician according to his operative needs
- Verify instrument integrity before reutilization: absence of small cracks, breakages and/or deformations on plastic part, diamond presence on metal band
- Be careful not to injure the soft tissues, especially when you cross the equator of the tooth
- To ensure traceability during the instrument lifetime we recommend to preserve the packaging
- On the product packages SDC lists the references that allow you to trace the entire production cycle of each product. In case of detected problems of a constructive nature of SDC product, we ask you to return in its full packaging, in order to allow us to be able to trace the cause of the problem that we will report

8. Maintenance

- The instrument is supplied non-sterile: sterilize before the first use
- In order to preserve the abrasive power of the file, after each use the instrument must be cleaned of any residues, disinfected and sterilized

9. Cleaning

- Disinfection can take place in an ultrasonic cleaning device containing a liquid degreaser and / or a chemically inhibited disinfectant
- The cleaning of the instrument after each use can be done by brushing the diamond part provided that the brush bristles are not too rigid to avoid damaging metal file, diamond and polymer

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10. Sterilization

- Use single or fractionated vacuum autoclave and vacuum drying
- The hot air sterilizers are not suitable for diamond instruments
- Sterilization must follow validated procedures in compliance with European standard EN 13060:
 - B Standard 134 (T=135.5°C, P=2.16bar, t_{Plateau}=4', t_{Drying}=15', t_{Total}=30'÷40')
 - B Prion 134 (T=135.5°C, P=2.16bar, t_{Plateau}=18', t_{Drying}=15', t_{Total}=44'÷54')
 - B Standard 121 (T=122.5°C, P=1.14bar, t_{Plateau}=15', t_{Drying}=20', t_{Total}=50'÷60')

11. Contraindications

- Use only as described above to avoid incurring drawbacks due to the carelessness and distraction
- Do not use on patients suffering from indisposition and / or diseases of any nature and / or type
- Do not disperse in the environment
- The instrument gets in contact with biological tissues during its usage: dispose of the instrument in accordance with the destination country legislation
- Keep out of the sight and reach of children
- Do not swallow

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